740 nm edge BrightLine® single-edge standard epi-fluorescence dichroic beamsplitter





Semrock, Inc

3625 Buffalo Road, Suite 6 Rochester, New York 14624

Main Phone: +1 585.594.7050 (worldwide)
Toll Free Phone: 866.736.7625 (866-SEMROCK)
(within US and Canada)

Your filter spectrum may differ slightly from the typical spectrum above, but is certified to meet the optical specifications noted below.



740 nm edge BrightLine® single-edge standard epi-fluorescence dichroic beamsplitter

Semrock offers a wide range of dichroic beamsplitters that exhibit steep edges with very high reflection and transmission bands. This standard epi-fluorescence dichroic beamsplitter is optimized only for broadband light sources such as LEDs and metal halide lamps.

Part Number	Size	Price1	Stock Status
FF740-Di01-25x36	25.2 mm x 35.6 mm x 1.1 mm (unmounted)	\$335	2nd Day Ship

Don't see a size you need? Contact us for custom sizing - delivery confirmed ARO (sizing fee applies).

1) US domestic pricing only. If you are ordering from outside the US, please contact your nearest regional distributor for the correct list price.

Optical Specifications

Specification	Value	
Reflection Band 1	Ravg > 98% 480 – 720 nm	
Edge Wavelength 1	740 nm	
Transmission Band 1	Tavg > 90% 750 - 825 nm	

General Filter Specifications

Specification	Value	
Angle of Incidence	45 ± 1.5 degrees	
Cone Half-angle	2 degrees	
Optical Damage Rating	Testing has proven to show no signs of degradation when exposed to at least 6.0 W of power from an unfiltered xenon arc lamp over a 25 mm diameter (corresponding to 1.2 W/cm²) for over 500 hrs.	
Flatness	Standard	
Steepness	Standard	
Filter Effective Index	1.62 Understanding 'Effective Index of Refraction' neff	
Flatness / RWE Classification	Standard Epi-fluorescence	
Reflected Wavefront Error	>> 6λ P-V RWE @ 632.8 nm	

Physical Filter Specifications (applies to standard sized parts; contact us regarding other sizes)

· · · · · · · · · · · · · · · · · · ·			
Specification	Value		
Transverse Dimensions (L x W)	25.2 mm x 35.6 mm		
Transverse Tolerance	± 0.1 mm		
Filter Thickness (unmounted)	1.05 mm		
Filter Thickness Tolerance (unmounted)	± 0.05 mm		
Clear Aperture	≥ 80% (elliptical)		
Scratch-Dig	60-40		

Substrate Thickness (unmounted)

Substrate Thickness Tolerance (unmounted)

± 0.05 mm

Orientation

Reflective surface marked with part number - Orient in direction of incoming light