
527/20 nm BrightLine® single-band bandpass filter

Part Number: FF01-527/20-25





Semrock, Inc.

3625 Buffalo Road, Suite 6 Rochester, New York 14624

Main Phone: +1 585.594.7050 (worldwide)
Toll Free Phone: 866.736.7625 (866-SEMROCK)
(within US and Canada)

Your filter spectrum may differ slightly from the typical spectrum above, but is certified to meet the optical specifications noted below.



527/20 nm BrightLine® single-band bandpass filter

Individual fluorescence bandpass filters that have been optimized for use in a variety of fluorescence instruments. All thin-film, hard-coated construction for unsurpassed performance and reliability.

Part Number	Size	Price1	Stock Status
FF01-527/20-25	25 mm x 3.5 mm	\$355	2 nd Day Ship
FF01-527/20-32	32 mm x 3.5 mm	\$582	2nd Day Ship
FF01-527/20-23.3-D	23.3 mm x 2.0 mm (unmounted)	\$355	2nd Day Ship

Don't see a size you need? Contact us for custom sizing - delivery confirmed ARO (sizing fee applies).

1) US domestic pricing only. If you are ordering from outside the US, please contact your nearest regional distributor for the correct list price.

Optical Specifications

Value	
Tavg > 93% 517 – 537 nm	
527 nm	
20 nm	
25.3 nm	
ODavg > 4 200 - 484 nm	
ODavg > 10 484 - 504 nm (Design specification - measurements are noise-floor limited)	
OD > 3.5 510.5 nm	
ODavg > 5 548 - 925 nm	
ODavg > 2 925 - 1100 nm	

General Filter Specifications

Specification	Value	
Angle of Incidence	0 ± 5 degrees	
Cone Half-angle	7 degrees	
Optical Damage Rating	Testing has proven to show no signs of degradation when exposed to at least 6.0 W of power from an unfiltered xenon arc lamp over a 25 mm diameter (corresponding to 1.2 W/cm²) for over 500 hrs.	
Filter Effective Index	1.86 Understanding 'Effective Index of Refraction' neff	

Physical Filter Specifications (applies to standard sized parts; contact us regarding other sizes)

Specification	Value
Transverse Dimensions (Diameter)	25 mm
Transverse Tolerance (mounted)	+ 0.0 / – 0.1 mm

Filter Thickness (Mounted)	3.5 mm
Filter Thickness Tolerance (Mounted)	± 0.1 mm
Clear Aperture	≥ 22 mm
Scratch-Dig	60-40
Substrate Thickness (unmounted)	2.0 mm
Substrate Thickness Tolerance (unmounted)	± 0.1 mm
Orientation Arrow on ring indicates preferred direction of propagation of light	