
#### 447/60 nm BrightLine® single-band bandpass filter

# Part Number: FF02-447/60-25





#### Semrock, Inc

3625 Buffalo Road, Suite 6 Rochester, New York 14624

Main Phone: +1 585.594.7050 (worldwide)
Toll Free Phone: 866.736.7625 (866-SEMROCK)
(within US and Canada)

Your filter spectrum may differ slightly from the typical spectrum above, but is certified to meet the optical specifications noted below.



### 447/60 nm BrightLine® single-band bandpass filter

Individual fluorescence bandpass filters that have been optimized for use in a variety of fluorescence instruments. All thin-film, hard-coated construction for unsurpassed performance and reliability.

Part Number	Size	Price1	Stock Status
FF02-447/60-25	25 mm x 3.5 mm	\$325	In Stock
FF02-447/60-25-STR	25 mm threaded ring for Sutter Lambda filter wheel	\$345	2nd Day Ship
FF02-447/60-32	32 mm x 3.5 mm	\$532	In Stock
FF02-447/60-23.3-D	23.3 mm x 2.0 mm (unmounted)	\$325	2nd Day Ship

Don't see a size you need? Contact us for custom sizing - delivery confirmed ARO (sizing fee applies).

1) US domestic pricing only. If you are ordering from outside the US, please contact your nearest regional distributor for the correct list price.

## Optical Specifications

Value
Tavg > 93% 417 – 477 nm
447 nm
60 nm
65.3 nm
ODavg > 2 200 - 352 nm
ODavg > 10 352 - 402 nm (Design specification - measurements are noise-floor limited)
OD > 4.5 409.5 nm
ODavg > 5 488 - 700 nm
ODavg > 5 700 - 925 nm
ODavg > 2 925 - 1100 nm

### **General Filter Specifications**

Specification	Value	
Angle of Incidence	0 ± 5 degrees	
Cone Half-angle	7 degrees	
Optical Damage Rating	Testing has proven to show no signs of degradation when exposed to at least 6.0 W of power from an unfiltered xenon arc lamp over a 25 mm diameter (corresponding to 1.2 W/cm²) for over 500 hrs.	
Filter Effective Index	2.18 Understanding 'Effective Index of Refraction' neff	

#### Physical Filter Specifications (applies to standard sized parts; contact us regarding other sizes)

Specification	Value
---------------	-------

Transverse Dimensions (Diameter)	25 mm
Transverse Tolerance (mounted)	+ 0.0 / - 0.1 mm
Filter Thickness (Mounted)	3.5 mm
Filter Thickness Tolerance (Mounted)	± 0.1 mm
Clear Aperture	≥ 22 mm
Scratch-Dig	60-40
Substrate Thickness (unmounted)	2.0 mm
Substrate Thickness Tolerance (unmounted)	± 0.1 mm
Orientation	Arrow on ring indicates preferred direction of propagation of light